Measurement of Density and Structural Short Range Order of Levitated Liquid Metals¹

G. Jacobs^{2,3}, I. Egry², E. Gorges², and M. Langen²

- Paper presented at the Thirteenth Symposium on Thermophysical Properties,

 June 22-27, 1997, Boulder, Colorado, U.S.A.
- Institut für Raumsimulation, Deutsche Forschungsanstalt für Luft- und Raumfahrt, D-51170 Köln, Germany.
- To whom correspondence should be addressed.

ABSTRACT

The determination of thermophysical properties and structure of undercooled metallic melts must be accomplished by a contactless method due to the high reactivity of the material. It has been shown that electromagnetic levitation provides high purity conditions to allow deep undercooling. Density and thermal expansion of a levitated drop can be derived from volume measurements using a CCD camera and a digital image processing system. Combining levitation with the extended x-ray absorption fine structure (EXAFS) spectroscopy leads to the possibility to study the local structure of the liquid in a wide temperature range including the deeply undercooled regime.

KEY WORDS: density; EXAFS; levitation; liquid metals; local structure; thermal expansion; undercooled melt.

1. INTRODUCTION

Information about thermophysical properties of liquid metals is important for materials science as well as for industrial processes and, in combination with knowledge about the structure of the melt, can help to understand the nature of the liquid state.

Due to their high temperature and reactivity, liquid metals are quite difficult to handle by conventional techniques. For this type of material, electromagnetic levitation is an elegant method, which allows containerless processing and, due to the high-purity environment, a deep undercooling of the melt. This non-equilibrium, metastable state is of interest for many reasons: First, a liquid is more quiescent in the undercooled region than it is in equilibrium, because the temperature is lower by several hundreds of degrees, and, consequently, thermal fluctuations are greatly reduced. This quiescence makes an undercooled liquid suitable for structural studies. In addition, an undercooled liquid is not in equilibrium and new, metastable, phases, normally excluded by the phase diagram, may form.

For studies on undercooled samples, electromagnetic levitation has to be complemented by non-contact diagnostic tools. With an advanced optical method using a video camera and a digital image processing system, we are able to measure several thermodynamic properties [1,2] e. g. the density and thermal expansion of the liquid. The short range structure can be investigated by combining levitation and the extended x-ray absorption fine structure spectroscopy (EXAFS). EXAFS is an established tool for the analysis of structural parameters like the coordination number, the nearest-neighbor distance and its variance.

2. EXPERIMENTAL

2.1 Electromagnetic levitation

Levitation of electrically conducting samples can be achieved by placing the sample into a high-frequency alternating inhomogeneous electromagnetic field, produced by a specially designed levitation coil. This field B induces eddy currents in the sample, which interact with the field. Levitation is caused by the Lorentz force F, which is, to lowest order in a multipole expansion, proportional the gradient of the square of the field

$$F \propto \nabla B^2 \tag{1}$$

Simultaneously the material is heated due to ohmic losses of the induced currents. The power P absorbed by the sample is proportional to the square of the field

$$P \propto B^2 \tag{2}$$

The levitated sample is positioned in a potential well generated by the electromagnetic field and performs oscillations about its equilibrium position, with a frequency in the order of 2-5 Hz. At the same time liquid samples display free surface oscillations, whose restoring force is the surface tension. Using samples with a typical mass of about 1 g, the frequency of these oscillations is in the order of 40-50 Hz.

A sketch of the levitation facility is shown in Fig. 1. The sample is processed inside of an UHV chamber in a levitation coil, which is connected to a 6 kW RF generator with a frequency of about 260 kHz. The temperature of the sample can be controlled by convective cooling using He/H₂ gas and measured with a pyrometer.

2.2. The density

There are a number of techniques to measure the density of liquid systems; however, non of them is a non-contact method and applicable to undercooled liquids. This can be achieved by using a videographic method. Taking a picture of the sample allows to fit its shape by a series expansion in Legendre polynomials [3]

$$r(\mathbf{u}) = \sum_{l} \varepsilon_{l} P^{l}(\mathbf{u}) \tag{3}$$

where $\,u=\cos\vartheta$. Provided that the liquid sample has a rotational symmetry, the volume is then given by

$$V = \frac{2\pi}{3} \int_{-1}^{1} r^3(u) du$$
 (4)

and, finally, the density is obtained from

$$\rho = \frac{M}{V} \tag{5}$$

where M is the mass of the sample.

For the measurement of the density an infrared and two polarization filters are put in front of the chamber as shown in Fig. 1a, and a CCD camera takes video sequences of the sample, which are stored on a video tape.

The filter arrangement is used to keep the intensity on the CCD camera constant: The sample's brightness is a function of temperature. By controlling the angle between two polarization filters, the intensity of the picture is kept constant. Since the polarization

filters work only in the visible range, we use an infrared filter with a very low transmission rate in the infrared wavelength regime.

The major experimental difficulties lie in the high spatial resolution required to resolve volume changes of the order of $\Delta V/V \approx 10^{-4}$ and in the elimination of the non-rotationally symmetric surface oscillations [4]. Using an image processing system and averaging over at least 100 pictures allows to evaluate the shape of the sample with sufficient accuracy.

2.3. EXAFS

The absorption coefficient of an atom in condensed matter is known to oscillate at energies of approximately 50-1000 eV above the absorption edge. This structure is called EXAFS and reflects the scattering of the emitted photoelectron by the surrounding atoms [5]. The EXAFS signal is defined as the normalized difference between the actual absorption $\alpha(k)$ and the absorption of an isolated atom $\alpha_0(k)$:

$$\chi(\mathbf{k}) = \frac{\alpha(\mathbf{k}) - \alpha_0(\mathbf{k})}{\alpha_0(\mathbf{k})} \tag{6}$$

where k is the wave vector of the emitted photoelectron.

In the case of a purely Gaussian pair distribution function the expression for $\chi(k)$ is given by [6]

$$\chi(k) = -\sum_{j} \frac{N_{j}}{R_{j}^{2} k} \left| f_{j}(\pi, k) \right| e^{-2\sigma_{j}^{2} k^{2}} e^{-2R_{j}/\lambda(k)} \sin[2kR_{j} + \phi_{j}(k)]$$
(7)

where N_j is the number of atoms in the j th shell, R_j is the mean distance between the absorber and the j th scatterer, σ_j^2 is the corresponding mean square displacement of the Gaussian pair correlation function, $\left|f_j(\pi,k)\right|$ is the absolute value of the scattering amplitude, $\varphi(k)$ is the total phase shift of the electron wave function, and $\lambda(k)$ is the electron mean free path.

For the EXAFS measurements the chamber is equipped with a Kapton window to allow the incident synchrotron beam to reach the sample (Fig. 1b). Since it is not possible to measure the absorption of the sample directly in transmission, due to the thickness of the sample (about 5 mm in diameter), we detect the secondary fluorescence by two photodiodes which are placed opposite to the direction of the incident beam.

Two serious problems affect the measurement as well as the data evaluation: First, the absorption of a levitating liquid, performing translational and surface oscillations, has to be detected with sufficient precision, and secondly, the thermal vibrations of the atoms together with the absence of a crystalline order lead to a huge damping of the EXAFS oscillations, due to the Debye-Waller factor in Eqn. 7, which makes it difficult to obtain sharply peaked spectra.

3. RESULTS

3.1. Density

Fig 2 shows the shape of a levitated droplet. The full line was produced by fitting a series expansion in Legendre polynomials up to the order of l=4. The agreement of the fit with the data points is obvious and allows to evaluate the density and thermal expansion of liquids in convincing accuracy.

The result of the density measurements on liquid $Cu_{20}Ni_{80}$ is shown in Fig 3. The solid line represents a least squares fit to the data points, whereas the dashed line shows the results of Watanabe [7]. It was possible to measure the density in a temperature range of more than 300°C, including an undercooling of about 230°C.

3.2. EXAFS

The first EXAFS measurements on levitated metallic melts were performed using liquid $Co_{80}Pd_{20}$ samples [8,9]. These measurements, performed at DESY, Hamburg, showed the possibility of obtaining reasonable EXAFS spectra of levitating melts [10].

Fig. 4 shows two EXAFS spectra of an overheated and an undercooled liquid sample obtained during recent measurements at the European Synchrotron Radiation Facility (ESRF, Grenoble/France) using a focused x-ray beam. By using the GNXAS method [11,12], it was possible to evaluate the mean distance of the Cobalt atoms to each other, $R_{\text{Co-Co}}$, as well as of the Cobalt atoms to the Palladium atoms, $R_{\text{Co-Pd}}$ (Fig. 5) and the corresponding Debye Waller factors $\sigma_{\text{Co-Co}}^2$ and $\sigma_{\text{Co-Pd}}^2$ (Fig. 6) in a temperature range of 420°C including an undercooling of more than 310°C.

3.3. Summary and Outlook

As we have shown, electromagnetic levitation opens the way to structural and thermophysical investigations of undercooled liquid metals. We plan to perform density measurements on CoPd. Combining the results of density measurements with the knowledge about the short range order in the liquid could lead to interesting insights into the origin of thermal expansion in liquids.

ACKNOWLEDGMENTS

Two of the authors (I. E. and G. J.) want to acknowledge A. Filipponi (BL 18, ESRF) for his help during the EXAFS measurements and the introduction to the GNXAS method.

REFERENCES

- I. Egry, S. Sauerland, G. Jacobs, *High Temperatures High Pressures* **26**: 217 (1994).
- I. Egry, G. Lohöfer, E. Gorges, G. Jacobs, *J. Phys.: Condens. Matter* **8**: 9363 (1996).
- 3 L. M. Racz, I. Egry, Rev. Sci. Instrum. 66: 4254 (1995).
- 4 E. Gorges, L. M. Racz, A. Schillings, I. Egry, *Int. J. Thermophys.* **17**: 1163 (1996).
- 5 B. K. Teo, in *EXAFS Spectroscopy*, edited by B. K. Teo and D. C. Joy (Plenum, New York, 1981), p. 13.
- P. Lagarde, in *Amorphous Solids and the Liquid State*, edited by N. H. March, R.A. Street, and M. Tosi (Plenum, New York, 1985), p. 365.
- 7 S. Watanabe, A. Saito, Nippon Kinzoku Gakkaishi **35**, 554 (1971)
- 8 D. Platzek, C. Notthoff, D. M. Herlach, G. Jacobs, D. Herlach, K. Maier, *Appl. Phys. Lett.* **65**: 1723 (1994).
- J. Reske, D. M. Herlach, F. Keuser, K. Maier, D. Platzek, *Phys. Rev. Lett.* 75, 737 (1995).
- 10 G. Jacobs, I. Egry, K. Maier, D. Platzek, J. Reske, R. Frahm, *Rev. Sci. Instrum*.67: 3683 (1996)
- 11 A. Filipponi, A. Di Cicco, C. R. Natoli, *Phys. Rev. B* **52**: 15122 (1995).
- 12 A. Filipponi, A. Di Cicco, *Phys. Rev. B* **52**: 15135 (1995).

- Fig. 1. Experimental setup for levitation experiments: a) optical system for density measurements; b) fluorescence radiation detection for EXAFS.
- Fig. 2. Shape of a levitated droplet: The solid line represents a fit to the data points.
- Fig. 3. Density of liquid $Cu_{20}Ni_{80}$. The solid line represents a least squares fit to the data points and the dashed line represents the literature value.
- Fig. 4. EXAFS spectra of liquid Co₈₀Pd₂₀ obtained at an overheating of about 110°C respectively an undercooling of about 290°C. The solid lines represent least squares fits to the data points. For reasons of clarity both curves are shifted against each other.
- Fig. 5. Mean value of the nearest neighbor distances $R_{\text{Co-Co}}$ and $R_{\text{Co-Pd}}$ in liquid $Co_{80}Pd_{20}$. The liquidus temperature is at 1337°C.
- Fig. 6. Debye Waller factors $\sigma_{\text{Co-Co}}$ and $\sigma_{\text{Co-Pd}}$ in liquid $\text{Co}_{80}\text{Pd}_{20}$. The liquidus temperature is at 1337°C.

Fig. 1a: Density













